



RAC-gamma serine, threonine-protein kinase(AKT3)

OVERVIEW

PRODUCT REFERENCE	PX-P4707
ORIGIN SPECIES	
FRAGMENT	Met1-Lys143
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	16.98 kDa
FORM	
BUFFER	PBS pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9Y243
ALIAS/SYNONYMS	Protein kinase Akt-3,Protein kinase B gamma,PKB gamma,RAC-PK-gamma,STK-2

