



Alpha-synuclein(SNCA)

OVERVIEW

PRODUCT REFERENCE	PX-P4717
ORIGIN SPECIES	
FRAGMENT	Met1-Ala140
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	14.46 kDa
FORM	
BUFFER	PBS pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P37840
ALIAS/SYNONYMS	Non-A beta component of AD amyloid,Non-A4 component of amyloid precursor,NACP,NACP, PARK1

