



## Receptor tyrosine-protein kinase erbB-3(ERBB3)Fc Tag

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P4725
<b>ORIGIN SPECIES</b>	
<b>FRAGMENT</b>	Met1-Thr643
<b>TAG</b>	C-terminal Fc Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Eukaryotic expression
<b>HOST</b>	Mammalian cells
<b>MOLECULAR WEIGHT</b>	97.56 kDa
<b>FORM</b>	
<b>BUFFER</b>	PBS pH 7.5
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	P21860
<b>ALIAS/SYNONYMS</b>	Proto-oncogene-like protein c-ErbB-3, Tyrosine kinase-type cell surface receptor HER3 , HER3

