



Pentraxin-related protein PTX3(PTX3)

OVERVIEW

PRODUCT REFERENCE	PX-P4726
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Ser381
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	41.98 kDa
FORM	
BUFFER	PBS pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P26022
ALIAS/SYNONYMS	Pentraxin-related protein PTX3, Tumor necrosis factor alpha-induced protein 5, TNF alpha-induced protein 5, Tumor necrosis factor-inducible gene 14 protein, TSG-14, TNFAIP5, TSG14

