



## Mucin-1(Ser24-Pro146)

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P4730
<b>ORIGIN SPECIES</b>	Homo sapiens (Human)
<b>FRAGMENT</b>	Ser24-Pro146
<b>TAG</b>	N-terminal His Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>HOST</b>	Escherichia coli (E.coli)
<b>MOLECULAR WEIGHT</b>	14.25 kDa
<b>FORM</b>	
<b>BUFFER</b>	PBS pH 7.5
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	P15941
<b>ALIAS/SYNONYMS</b>	MUC-1,Breast carcinoma-associated antigen DF3,Cancer antigen 15-3,CA 15-3,Carcinoma-associated mucin,Episialin,H23AG,Krebs von den Lungen-6,KL-6 ,PEMT,Peanut-reactive urinary mucin,PUM,Polymorphic epithelial mucin,PEM,Tumor-associated epithelial membrane

