



KAT8 regulatory NSL complex subunit 1(Kansl1)

OVERVIEW

PRODUCT REFERENCE	PX-P4737
ORIGIN SPECIES	
FRAGMENT	Lys814-Arg1036
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	25.60kDa
FORM	
BUFFER	PBS pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q80TG1
ALIAS/SYNONYMS	NSL complex protein NSL1, Non-specific lethal 1 homolog, Kiaa1267, Nsl1

