



c-JUN Protein - Transcription factor AP-1(JUN)

OVERVIEW

PRODUCT REFERENCE	PX-P4739
ORIGIN SPECIES	
FRAGMENT	Met1-Phe331
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	37.78kDa
FORM	
BUFFER	0.02% Sarcosyl, PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P05412
ALIAS/SYNONYMS	Activator protein 1,Proto-oncogene c-Jun,V-jun avian sarcoma virus 17 oncogene homolog,p39

