

Medium-chain specific acyl-CoA dehydrogenase, mitochondrial(Acadm)

OVERVIEW

PRODUCT REFERENCE	PX-P4740
ORIGIN SPECIES	
FRAGMENT	Lys26-Asn421
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	44.78kDa
FORM	
BUFFER	PBS pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P45952
ALIAS/SYNONYMS	MCAD, Acadm