

Electron transfer flavoprotein subunit alpha, mitochondrial(Etfa)

OVERVIEW

PRODUCT REFERENCE	PX-P4741
ORIGIN SPECIES	
FRAGMENT	Gln20-Lys333
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	33.99kDa
FORM	
BUFFER	50mM Tris-HCl pH 8.0,150mMNaCl
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q99LC5
ALIAS/SYNONYMS	Alpha-ETF, Etfa