

Cytochrome P450 1B1(CYP1B1)GST Tag

OVERVIEW

PRODUCT REFERENCE	PX-P4746
ORIGIN SPECIES	
FRAGMENT	Leu224-Asp492
TAG	N-Terminal GST Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	57.03kDa
FORM	
BUFFER	TBS, pH 8.0, 0.02%Sarcosyl
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q16678
ALIAS/SYNONYMS	CYP1B1,Hydroperoxy icosatetraenoate dehydratase