



Unconventional prefoldin RPB5 interactor 1(His Tag)

OVERVIEW

| | |
|---------------------------|--------------------|
| PRODUCT REFERENCE | PX-P4749 |
| ORIGIN SPECIES | |
| FRAGMENT | Met1-Asp535 |
| TAG | C-terminal His Tag |
| DELIVERY CONDITION | Dry Ice |

PRODUCT INFORMATION

| | |
|--------------------------------|--|
| EXPRESSION SYSTEM | Eukaryotic expression |
| HOST | Mammalian cells |
| MOLECULAR WEIGHT | 62.92kDa |
| FORM | |
| BUFFER | PBS, pH 7.5. |
| STABILITY & STORAGE | 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection) |
| UNIPROT | O94763 |
| ALIAS/SYNONYMS | C19orf2, NNX3, PPP1R19, RMP, URI |

