

10 kDa culture filtrate protein (esxB)GST Tag

OVERVIEW

PRODUCT REFERENCE	PX-P4754
ORIGIN SPECIES	
FRAGMENT	Met1-Phe100
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	37.27kDa
FORM	
BUFFER	TBS, pH8.0
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P9WNK4
ALIAS/SYNONYMS	ESXB_MYCTO, ESAT-6-like protein EsxB, CFP-10, Secreted antigenic protein MTSA-10