



RELA Protein - Transcription factor p65(RELA)

OVERVIEW

PRODUCT REFERENCE	PX-P4763
ORIGIN SPECIES	
FRAGMENT	Pro19-Asp291
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	32.22kDa
FORM	
BUFFER	PBS PH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q04206
ALIAS/SYNONYMS	Nuclear factor NF-kappa-B p65 subunit,Nuclear factor of kappa light polypeptide gene enhancer in B-cells 3,NFKB3

