

Mucin-4(MUC4)(2078-2117)

OVERVIEW

PRODUCT REFERENCE	PX-P4770
ORIGIN SPECIES	
FRAGMENT	Val2078-Glu2117
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	30.58kDa
FORM	
BUFFER	TBS, pH8.0
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q99102
ALIAS/SYNONYMS	MUC-4, Ascites sialoglycoprotein, ASGP, Pancreatic adenocarcinoma mucin, Testis mucin, Tracheobronchial mucin

