



## Thioredoxin domain-containing protein 5(TXNDC5)

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P4776
<b>ORIGIN SPECIES</b>	
<b>FRAGMENT</b>	Met1-Leu432
<b>TAG</b>	
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Eukaryotic expression
<b>HOST</b>	Mammalian cells
<b>MOLECULAR WEIGHT</b>	47,55
<b>FORM</b>	
<b>BUFFER</b>	PBS pH 7.5
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	Q8NBS9
<b>ALIAS/SYNONYMS</b>	EndoPDI; ERp46; FLJ90810; HCC-2; member 15; MGC3178; PDIA15; thioredoxin domain containing 5 (endoplasmic reticulum); thioredoxin domain containing 5; thioredoxin domain-containing protein 5; thioredoxin related protein; TLP46; TXNDC5

