



Vascular cell adhesion protein 1(VCAM1)

OVERVIEW

PRODUCT REFERENCE	PX-P4786
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Phe25-Glu698
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	74.24 kDa
FORM	
BUFFER	PBS pH 7.5, 0.02% SKL
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P19320
ALIAS/SYNONYMS	V-CAM 1, VCAM-1, INCAM-100, CD106

