



Endoplasmin(HSP90B1)

OVERVIEW

PRODUCT REFERENCE	PX-P4788
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Asp22-Ala311
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	32.62 kDa
FORM	
BUFFER	PBS pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P14625
ALIAS/SYNONYMS	94 kDa glucose-regulated protein,GRP-94,Heat shock protein 90 kDa beta member 1,Tumor rejection antigen 1,gp96 homolog,GRP94, TRA1

