



Annexin A2(ANXA2)

OVERVIEW

PRODUCT REFERENCE	PX-P4794
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Asp339
TAG	N terminus His tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	38.61 kDa
FORM	
BUFFER	PBS pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P07355
ALIAS/SYNONYMS	Annexin II, Annexin-2, Calpactin I heavy chain, Calpactin-1 heavy chain, Chromobindin-8, Lipocortin II, Placental anticoagulant protein IV, PAP-IV, Protein I, p36, ANX2, ANX2L4, CAL1H, LPC2D

