



## Estrogen receptor (in E.coli)

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P4797
<b>ORIGIN SPECIES</b>	Homo sapiens (Human)
<b>FRAGMENT</b>	Met297-Ser554
<b>TAG</b>	N terminus His tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>HOST</b>	Escherichia coli (E.coli)
<b>MOLECULAR WEIGHT</b>	29.47 kDa
<b>FORM</b>	
<b>BUFFER</b>	PBS pH 7.5, 0.02% SKL
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	P03372
<b>ALIAS/SYNONYMS</b>	ER, ER-alpha, Estradiol receptor, Nuclear receptor subfamily 3 group A member 1, ESR, NR3A1

