

Integrin beta-1(ITGB1)

OVERVIEW

PRODUCT REFERENCE	PX-P4798
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Gln21-Asp728
TAG	N terminus His tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	78.36 kDa
FORM	
BUFFER	PBS pH 7.5, 0.02% SKL
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P05556
ALIAS/SYNONYMS	Fibronectin receptor subunit beta, Glycoprotein IIa, GPIIA, VLA-4 subunit beta, CD_antigen: CD29, FNRB, MDF2, MSK12

