

60S ribosomal protein L10-like(RPL10L)

OVERVIEW

PRODUCT REFERENCE	PX-P4800
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Ser214
TAG	N terminus His tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	24.52 kDa
FORM	
BUFFER	PBS pH 7.5, 0.02% SKL
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P27635
ALIAS/SYNONYMS	DXS648E, QM,Laminin receptor homolog, Large ribosomal subunit protein uL16, Protein QM, Ribosomal protein L10, Tumor suppressor QM

