



Cytochrome c oxidase assembly factor 8(COA8)

OVERVIEW

PRODUCT REFERENCE	PX-P4817
ORIGIN SPECIES	
FRAGMENT	Ala40-Asn206
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	21.3kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q96IL0
ALIAS/SYNONYMS	APOP1, APOPT1, APOP-1, Apoptogenic protein 1, mitochondrial

