



# Proteogenix

Antibodies and Proteins

**Proteogenix SAS**  
Espace Européen de l'Entreprise  
19 rue de La Haye  
67300 Schiltigheim, France  
Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88  
contact@proteogenix.fr  
www.Proteogenix.science

## Cytochrome c oxidase assembly factor 8(COA8)

### OVERVIEW

---

<b>CATALOG NO</b>	PX-P4817
<b>ORIGIN SPECIES</b>	
<b>FRAGMENT</b>	Ala40-Asn206
<b>MOLECULAR WEIGHT</b>	21.3kDa
<b>DELIVERED WITH TAG</b>	N-terminal His Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>HOST</b>	Escherichia coli (E.coli)
<b>PURITY</b>	>90% by SDS-PAGE
<b>UNIPROT ID</b>	Q96IL0
<b>FORM</b>	Liquid
<b>BUFFER</b>	PBS, pH7.5
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

### More info

#### SYNONYMS

APOP1, APOPT1, APOP-1, Apoptogenic protein 1, mitochondrial

#### UNIPROT LINK

<https://www.uniprot.org/uniprot/Q96IL0>

### SEQUENCE

MGSHHHHHHSGAPERGAERRDTAPSGVSRFCPPRKSCHDWIGPPDKYSNLRPVHFIYIPENESPLEQKLRKLRQETQEWNNQQFWANQNLTFSKEKEEFIHSRLKTKGLGLRTESGQKATLNAEEMADFYKEFLSKNFQKHMYYNRDQWYKRNFAITFFMGKVALERIWNKLLKQKQKRSN