



## Interleukin-17A(IL17A)

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P4824
<b>ORIGIN SPECIES</b>	
<b>FRAGMENT</b>	Met1-Ala155
<b>TAG</b>	C-terminal His Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Eukaryotic expression
<b>HOST</b>	Mammalian cells
<b>MOLECULAR WEIGHT</b>	18.5kDa
<b>FORM</b>	
<b>BUFFER</b>	PBS, pH7.5
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	Q16552
<b>ALIAS/SYNONYMS</b>	Cytotoxic T-lymphocyte-associated antigen 8,CTLA-8,CTLA8, IL17,IL-17, IL-17A

