



Receptor-type tyrosine-protein phosphatase-like N(PTPRN)

OVERVIEW

PRODUCT REFERENCE	PX-P4827
ORIGIN SPECIES	
FRAGMENT	His603-Gln979
TAG	N-terminal MAT Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	44.8kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q16849
ALIAS/SYNONYMS	R-PTP-N, Islet cell antigen 512, ICA 512, Islet cell autoantigen 3, PTP IA-2, ICA3, ICA512

