



Putative trypsinogen(try)

OVERVIEW

PRODUCT REFERENCE	PX-P4830
ORIGIN SPECIES	
FRAGMENT	Asp20-Ser247
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Yeast
MOLECULAR WEIGHT	25.3kDa
FORM	
BUFFER	PBS, pH 6.0
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	C6L245
ALIAS/SYNONYMS	PRSS2

