



RNA-binding protein 38(RBM38)

OVERVIEW

PRODUCT REFERENCE	PX-P4832
ORIGIN SPECIES	
FRAGMENT	Thr132-Gln239
TAG	N-terminal GST Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	38.3kDa
FORM	
BUFFER	PBS, pH 7.5, 4M urea
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9H0Z9
ALIAS/SYNONYMS	CLL-associated antigen KW-5,HSRNASEB, RNA-binding motif protein 38, RNA-binding region-containing protein 1, ssDNA-binding protein SEB4, RNPC1, SEB4

