



NF-kappa-B inhibitor alpha(NFKBIA) (with CPP2)

OVERVIEW

PRODUCT REFERENCE	PX-P4835
ORIGIN SPECIES	
FRAGMENT	Met1-Leu317
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	40.1kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P25963
ALIAS/SYNONYMS	I-kappa-B-alpha, IκB-alpha, IκappaBα, Major histocompatibility complex enhancer-binding protein MAD3, IKBA, MAD3, NFKBI

