



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

T-cell activation Rho GTPase-activating protein(TAGAP)

OVERVIEW

PRODUCT REFERENCE	PX-P4838
ORIGIN SPECIES	
FRAGMENT	Met1-Ile731
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	81.8kDa
FORM	
BUFFER	PBS,pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q8N103
ALIAS/SYNONYMS	T-cell activation GTPase-activating protein,TAGAP1



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.