



## Indoleamine 2,3-dioxygenase 2(IDO2)

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P4851
<b>ORIGIN SPECIES</b>	
<b>FRAGMENT</b>	Met14-Gly420
<b>TAG</b>	N-terminal His Tag & sumo Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>HOST</b>	Escherichia coli (E.coli)
<b>MOLECULAR WEIGHT</b>	58.2kDa
<b>FORM</b>	
<b>BUFFER</b>	PBS, pH7.5
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	Q6ZQW0
<b>ALIAS/SYNONYMS</b>	Indoleamine 2,3-dioxygenase-like protein 1,Indoleamine-pyrrole 2,3-dioxygenase-like protein 1,IDO-2,INDOL1

