



## Pyruvate kinase PKM(PKM)

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P4858
<b>ORIGIN SPECIES</b>	
<b>FRAGMENT</b>	Met1-Pro531
<b>TAG</b>	N-terminal His Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>HOST</b>	Escherichia coli (E.coli)
<b>MOLECULAR WEIGHT</b>	60.1kDa
<b>FORM</b>	
<b>BUFFER</b>	PBS, pH7.5
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	P14618
<b>ALIAS/SYNONYMS</b>	Cytosolic thyroid hormone-binding protein,CTHBP,Opa-interacting protein 3,OIP-3,Pyruvate kinase 2/3,Pyruvate kinase muscle isozyme,Thyroid hormone-binding protein 1,THBP1,Tumor M2-PK,p58,OIP3, PK2, PK3, PKM2

