



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Receptor-interacting serine, threonine-protein kinase 4(RIPK4)

OVERVIEW

PRODUCT REFERENCE	PX-P4863
ORIGIN SPECIES	
FRAGMENT	Phe22-Ile286
TAG	N-terminal GST Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	57.6kDa
FORM	
BUFFER	
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P57078
ALIAS/SYNONYMS	ANKRD3, DIK, Ankyrin repeat domain-containing protein 3, PKC-delta-interacting protein kinase



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.