



Phosphatidylethanolamine-binding protein 1(PEBP1)

OVERVIEW

PRODUCT REFERENCE	PX-P4869
ORIGIN SPECIES	
FRAGMENT	Met1-Lys187
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	21.98kDa
FORM	
BUFFER	
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P30086
ALIAS/SYNONYMS	PEBP-1, HCNPpp, RKIP, HCNP, PBP, PEBP

