



Red fluorescent protein drFP583

OVERVIEW

PRODUCT REFERENCE	PX-P4872
ORIGIN SPECIES	
FRAGMENT	Met1-Leu225
TAG	N-terminal HN Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	30.98kDa
FORM	
BUFFER	
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9U6Y8
ALIAS/SYNONYMS	DsRed

