



Sodium-dependent phosphate transport protein 2B(SLC34A2)

OVERVIEW

PRODUCT REFERENCE	PX-P4875
ORIGIN SPECIES	
FRAGMENT	Gly191-Leu362
TAG	N-terminal GST Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	46.27kDa
FORM	
BUFFER	
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	O95436
ALIAS/SYNONYMS	Sodium-phosphate transport protein 2B,Na(+)-dependent phosphate cotransporter 2B,NaPi3b,Sodium/phosphate cotransporter 2B,Na(+)/Pi cotransporter 2B,NaPi-2b,Solute carrier family 34 member 2

