



Interleukin-34(IL34)

OVERVIEW

PRODUCT REFERENCE	PX-P4881
ORIGIN SPECIES	House Mouse
FRAGMENT	Asn21-Pro235
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	27.62kDa
FORM	
BUFFER	
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q8R1R4
ALIAS/SYNONYMS	IL-34, IL34

