



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Pro-neuregulin-1, membrane-bound isoform(NRG1)

OVERVIEW

PRODUCT REFERENCE	PX-P4883
ORIGIN SPECIES	
FRAGMENT	Lys22-Lys246
TAG	C-terminal Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	49.43kDa
FORM	
BUFFER	
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q02297-6
ALIAS/SYNONYMS	ARIA; GGF; GGF2; HGL; HRG; HRG1; HRGA; MST131; MSTP131; NDF; neuregulin 1; Neuregulin1 beta 1; Neuregulin-1 beta 1; NRG1-IT2; SMDF



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.