

Fibroblast Growth Factor proteins 1(FGF1) (16-155)

OVERVIEW

PRODUCT REFERENCE	PX-P4884
ORIGIN SPECIES	
FRAGMENT	Phe16~Asp155
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	17kDa
FORM	
BUFFER	
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P05230
ALIAS/SYNONYMS	FGFA,aFGF,ECGF,HBGF-1