

Human FAP protein

OVERVIEW

PRODUCT REFERENCE	PX-P4887
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Ile523-Asp760
TAG	His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	27.55 kDa
FORM	Frozen Liquid or Lyophilized
BUFFER	PBS,pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles)
UNIPROT	
ALIAS/SYNONYMS	Prolyl endopeptidase FAP,Prolyl endopeptidase FAP,Dipeptidyl peptidase FAP,Fibroblast activation protein alpha,FAPalpha,Gelatine degradation protease FAP,Integral membrane serine protease,Post-proline cleaving enzyme,Serine integral membrane protease,SIMP