



IdeS Protease

OVERVIEW

PRODUCT REFERENCE	PX-P4888
ORIGIN SPECIES	Streptococcus pyogenes
FRAGMENT	
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	
FORM	liquid
BUFFER	PBS pH 7.5, 50%glycerol
STABILITY & STORAGE	Store at -20 °C for 12 months from the date of receipt.
UNIPROT	
ALIAS/SYNONYMS	Immunoglobulin G-degrading enzyme of Streptococcus pyogenes

