



SUMO Protease

OVERVIEW

PRODUCT REFERENCE	PX-P4890
ORIGIN SPECIES	
FRAGMENT	
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	26.55kDa
FORM	Frozen Liquid or Lyophilized
BUFFER	50 mM Tris-HCl, pH 8.0, 0.2% Igepal (NP-40), 150 mM NaCl, 1mM DTT, 50% glycerol
STABILITY & STORAGE	Store at -20 ° C, 12 months shelf life
UNIPROT	
ALIAS/SYNONYMS	

