



# Proteogenix

**Proteogenix SAS**

www.proteogenix.science  
19 rue de la Haye  
67300 Schiltigheim, France  
Tel +33 (0) 390 20 54 70  
Fax +33 (0) 978 53 36 90  
products@proteogenix.fr

## Taq DNA Polymerase

### OVERVIEW

---

**PRODUCT REFERENCE** PX-P4892  
**ORIGIN SPECIES**  
**FRAGMENT**  
**TAG**  
**DELIVERY CONDITION** Dry Ice

### PRODUCT INFORMATION

---

**EXPRESSION SYSTEM** Prokaryotic expression  
**HOST** Escherichia coli (E.coli)  
**MOLECULAR WEIGHT** 94kDa  
**FORM** Frozen Liquid or Lyophilized  
**BUFFER** 20mM Tris-HCl (pH8.0), 0.1 mM EDTA, 1mM DTT, 100 mM KCl, Stabilizers, 50% glycerol  
**STABILITY & STORAGE** Store at -20 ° C, 12 months shelf life  
**UNIPROT**  
**ALIAS/SYNONYMS**



**WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.**