



SIRP alpha human recombinant protein (His-Tag)

OVERVIEW

PRODUCT REFERENCE	PX-P4897
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Lys32-Asn372
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	37.51kDa
FORM	
BUFFER	PBS,pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P78324
ALIAS/SYNONYMS	BIT, MFR, MYD1, PTPNS1, SHPS2, SIRP

