

Cow PAG1B Recombinant Protein (Partial)

OVERVIEW

PRODUCT REFERENCE	PX-P4898
ORIGIN SPECIES	Bos taurus
FRAGMENT	Partial
TAG	C-Ter His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	38,07kDa
FORM	liquid
BUFFER	PBS pH 7.5, 0.02% SKL
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q29432
ALIAS/SYNONYMS	

