



## Apolipoprotein ApoA1 recombinant protein, C-His Tag

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P4905
<b>ORIGIN SPECIES</b>	Homo sapiens (Human)
<b>FRAGMENT</b>	Arg19-Gln267
<b>TAG</b>	C-terminal His Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Eukaryotic expression
<b>HOST</b>	Mammalian cells
<b>MOLECULAR WEIGHT</b>	27,39kDa
<b>FORM</b>	
<b>BUFFER</b>	PBS,pH7.5
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	P02647
<b>ALIAS/SYNONYMS</b>	ApoA-1, ApoA1

