

Human BDNF

OVERVIEW

PRODUCT REFERENCE	PX-P4907
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Arg128-Arg247
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	15.16kDa
FORM	
BUFFER	PBS,pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P23560
ALIAS/SYNONYMS	Brain-derived neurotrophic factor,BDNF,Abrineurin,BDNF precursor form,ProBDNF

