



Stanniocalcin 1

OVERVIEW

PRODUCT REFERENCE	PX-P4987
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Ser28-Ala247
TAG	GST Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	49.49 kDa
FORM	
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P52823
ALIAS/SYNONYMS	STC1

