



LTA recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5015
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Leu35-Leu205
TAG	His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	20.57 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P01374
ALIAS/SYNONYMS	Lymphotoxin-alpha,LT-alpha,TNF-beta,Tumor necrosis factor ligand superfamily member 1, TNFb

