



Human C5 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5117
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Asp822-Lys925
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	188.33 kDa
FORM	
BUFFER	PBS,pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P01031
ALIAS/SYNONYMS	Complement C5, C3 and PZP-like alpha-2-macroglobulin domain-containing protein 4, C5 anaphylatoxin, Complement C5 alpha, C5 beta

