



Human DLL4 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5121
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Leu16-Ala539
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	59.27 kDa
FORM	
BUFFER	PBS,pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9NR61
ALIAS/SYNONYMS	Delta-like protein 4, Drosophila Delta homolog 4, Delta4, DLL4

